



Instrument Engineering Change Notice

NO. Glas-2102


Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek Glas DATE : December 8, 2021
PRODUCT CODES : N/A ORIGINATOR : Hirabayashi
IECN TYPE : NEW PART REFERENCE NUMBER : 20311
SERIAL NUMBERS AFFECTED : N/A ATTACHMENTS : N/A
FIELD RETROFIT REQUIRED : NO

Title: Addition of Alternative 20-Slide Basket

[REASON AND DESCRIPTION FOR CHANGE]

Due to difficulty in procuring base material for the current 20-slide baskets (4768), the manufacturer decided to use the alternative material (RENATUS®). With this change, the basket color is changed to white.



Left: alternative basket
Right: current basket

NOTE:

The quality and performance of alternative 20-slide baskets were proved equal to the current 20-slide baskets through evaluation testing by the manufacturer. The alternative basket is planned to be supplied from the beginning of the year 2022.

However, note that the change is just a temporary expedient until the current base material can be obtained sufficiently and stably.

[PART CHANGE] (IVDP=Inventory Disposal, ICHB=Inter Changeable)

| Description | Old | QTY | IVDP | New | QTY | ICHB |
|----------------------------------|-----|-----|------|------------------|-----|------|
| 20-Slide Basket (alternative) | --- | --- | | A-MK25-0495-00-A | --- | |

[REMARKS]

N/A

[DISTRIBUTION]

- SFA
- SFJ
- SFE
- SCN

[APPROVAL]

Name : T. Sugawara

Signature :

