



Instrument Engineering Change Notice

NO. Glas g2-2002

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek Glas g2 DATE : August 31, 2020
 PRODUCT CODES : 6500, 6501, 6502, 6499 ORIGINATOR : T. Iijima
 IECN TYPE : PART DESIGN CHANGE REFERENCE NUMBER : 16990
 SERIAL NUMBERS AFFECTED : See Remarks ATTACHMENTS : N/A
 FIELD RETROFIT REQUIRED : NO

Title: Hardware Change for New CI Logo and Color

[REASON AND DESCRIPTION FOR CHANGE]

The control panel and interior label are to be changed to the new CI logo and color.

[PART CHANGE] (IVDP=Inventory Disposal, ICHB=Inter Changeable)

Description	Old	QTY	IVDP	New	QTY	ICHB
Control Panel (Glas g2 CI) / 010-S19054	---	---		00001-4180	1	NO
Control Panel / 010-S07104(02)	A4150189	1	NO	---	---	NO
Interior Label (CI)	---	---		A-MN94-0770-00-A	1	NO
Interior Label	N94-762-01	1	NO	---	---	NO
Interior Label (CHN) (CI)	---	---		A-MN94-0771-00-A	1	NO
Interior Label (CHN)	A-MN94-0767-00-A	1	NO	---	---	NO

[REMARKS]

Implementation:

A1: 6500 0581 and after

J0: 6501 0731 and after

E2: 6502 0577 and after

JC2: 6499 0072 and after

[DISTRIBUTION]

SFA
 SFJ
 SFE
 SCN

[APPROVAL]

Name : T. Sugawara

Signature :

<REFERENCE>

Current Logo and Color

New CI Logo and Color

* Control Panel



* Interior Label

