



# Instrument Engineering Change Notice

NO. Glas g2-2003

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek Glas g2                      DATE : October 19, 2020  
PRODUCT CODES : 6500, 6501, 6502, 6499              ORIGINATOR : Takamura  
IECN TYPE : Part Design Change                      REFERENCE NUMBER : 16857  
SERIAL NUMBERS AFFECTED : See Remarks              ATTACHMENTS : N/A  
FIELD RETROFIT REQUIRED : NO

## Title: Change of Oil Damper

### [ REASON AND DESCRIPTION FOR CHANGE ]

Due to product discontinuation, the oil damper for opening and closing the loading door is to be replaced with the new one.

Because the new oil damper is not interchangeable with the old damper and needs to be installed in a different way, the parts to install the oil damper are also redesigned.

Therefore, the new oil damper needs to be replaced with the old one together with other new parts as a set.

### [ PART CHANGE ] ( IVDP=Inventory Disposal, ICHB=Inter Changeable )

Description	Old	QTY	IVDP	New	QTY	ICHB
Oil Damper B-466-3	D9280283	1	NO	---	---	NO
Oil Damper B-499-3	---	---		00001-3531	1	NO
Loading Door Damper Mounting Plate	N94-705-00	1 (*)	NO	---	---	NO
Loading Door Damper Mounting Plate A	---	---		A-MN94-0058-00-A	1	NO
Flange Bearing / JFM-0608-06	---	---		00001-3484	2	NO
Plastic Washer / PPSW6	---	---		00001-3532	2	NO

(\*) One of the two damper mounting plates currently used is replaced with the new damper mounting plate A.

### [ REMARKS ]

Implementation:

A1: 6500 0581 and after  
J0: 6501 0731 and after  
E2: 6502 0577 and after  
JC2: 6499 0072 and after

---

**[ DISTRIBUTION ]**

SFA

SFJ


SFE

SCN

**[ APPROVAL ]**

Name : T. Sugawara

Signature :

A handwritten signature in black ink, appearing to read 'T. Sugawara', is written over a vertical line. The signature is cursive and stylized.