



# Instrument Engineering Change Notice

NO. Glas g2-2104

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek Glas g2                      DATE : November 30, 2021  
 PRODUCT CODES : 6500, 6501, 6502, 6499              ORIGINATOR : Y. Yazawa  
 IECN TYPE : NEW PART                                      REFERENCE NUMBER : 17164-1  
 SERIAL NUMBERS AFFECTED : See Remarks              ATTACHMENTS : N/A  
 FIELD RETROFIT REQUIRED : NO

## Title: Change of Serial Number Labels (New CI)

### [ REASON AND DESCRIPTION FOR CHANGE ]

The serial number labels were modified to uniform label layout and font type for new corporate identity. Also, the ZIP code for SFJ was changed in some labels.

For details, see the figures of the new labels on the next page.

(NOTE) This IECN only notifies of change of the serial number labels. They are not supplied as a part for service/maintenance.

### [ PART CHANGE ] ( IVDP=Inventory Disposal, ICHB=Inter Changeable )

Description	Old	QTY	IVDP	New	QTY	ICHB
Serial Number Label (Glas g2-J0)	N94-770-00	1		N94-770-00	1	YES
Serial Number Label (Glas g2-A1)	N94-771-01	1		N94-771-02	1	YES
Serial Number Label (Glas g2-E2)	N94-772-01	1		N94-772-02	1	YES
Serial Number Label (Glas g2-JC2)	A-MN74-0777-00-A	1		A-MN94-0777-01-A	1	YES

### [ REMARKS ]

Implementation:

A1: 6500 0602 and after

J0: 6501 0772 and after

E2: 6502 0605 and after

JC2: 6499 0077 and after

### [ DISTRIBUTION ]

SFA   
 SFJ   
 SFE   
 SCN

### [ APPROVAL ]

Name : T. Sugawara

Signature :

## New Serial Number Labels

- A1 model

**Tissue-Tek® Glas™ g2**  
Automated Glass Coverslipper

CODE	MODEL	VOLTS	AMPS	HERTZ
6500	Glas g2-A1	115	2.2	50/60

Serial Number : 65000000-0000

Manufactured for:  
Sakura Finetek USA, Inc., Torrance, CA 90501 U.S.A.  
Sakura Finetek Europe B.V., 2408 AV Alphen aan den Rijn, NL.  
Sakura Finetek Japan Co., Ltd., Tokyo, 103-0007, Japan

Manufactured by:  
Sakura Seiki Co., Ltd., Nagano, 387-0015, Japan

N94-771-02  
**MADE IN JAPAN**

**IVD**

- J0 model

**ティッシュ・テック® グラス™ ジー2** 6501

型 式	Glas g2-J0		
電 源 入 力	2.5A/AC100V 50/60Hz		
製造販売届出番号	20B2X00014000023	一般医療機器 (クラスI)	
製 造 番 号	65010000-0000	特定保守管理医療機器	

製造販売業者 **サクラ精機株式会社**  
長野県千曲市大字八幡 1122-8

販売業者 **サクラファインテックジャパン株式会社**  
東京都中央区日本橋本町 3-1-9

N94-770-01

**EMC適合**

- E2 model

**Tissue-Tek® Glas™ g2**  
Automated Glass Coverslipper

CODE	MODEL	VOLTS	AMPS	HERTZ
6502	Glas g2-E2	230	1.2	50/60

Serial Number : 65020000-0000

Manufactured for:  
Sakura Finetek USA, Inc., Torrance, CA 90501 U.S.A.  
Sakura Finetek Europe B.V., 2408 AV Alphen aan den Rijn, NL.  
Sakura Finetek Japan Co., Ltd., Tokyo, 103-0007, Japan

Manufactured by:  
Sakura Seiki Co., Ltd., Nagano, 387-0015, Japan

N94-772-02  
**MADE IN JAPAN**

**IVD**  
**CE**

- JC2 model

**自动盖片机**

<ul style="list-style-type: none"> <li>■ 产品名称: 自动盖片机</li> <li>■ 产品编号: 64990000-0000</li> <li>■ 规格型号: Glas g2-JC2</li> <li>■ 尺寸: 750(W)×620(D)×750(H) mm</li> <li>■ 生产者名称: Sakura Seiki Co., Ltd.</li> <li>■ 生产日期: 0000年 00月</li> <li>■ 有效期: 八年</li> </ul>	<ul style="list-style-type: none"> <li>■ 生产者地址: 75-5, IMOJIYA, CHIKUMA-CITY, NAGANO PREF., JAPAN</li> <li>■ 生产场所地址: 75-5, IMOJIYA, CHIKUMA-CITY, NAGANO PREF., JAPAN</li> <li>■ 电源连接条件: <ul style="list-style-type: none"> <li>● 电流: 5A</li> <li>● 电压: 230VAC ± 10%, 50/60Hz</li> <li>● 湿度: 30-85% (无冷凝)</li> <li>● 温度: 10-40℃</li> </ul> </li> </ul>
---	---

※其他内容详见说明书

 代理人: 樱花医疗科技(泰州)有限公司  
地 址: 泰州市药城大道1号R09幢  
邮编: 225300 电话: 0523-86201011 传真: 0523-86201022

MN94-0777-01