



# Instrument Engineering Change Notice

NO. TEC6-2202

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME	: Tissue-Tek TEC6	DATE	: June 15, 2022
PRODUCT CODES	: 5104, 5108, 5112, 5116, 5128, 5105, 5109, 5113, 5117, 5129	ORIGINATOR	: Yoshizawa
IECN TYPE	: PART DESIGN CHANGE	REFERENCE NUMBER	: 20017, 20020, 20122-1, 20123-1
SERIAL NUMBERS AFFECTED	: See Remarks	ATTACHMENTS	: N/A
FIELD RETROFIT REQUIRED	: NO		

## Title: Appearance Change of Circuit Protector Units

### [ REASON AND DESCRIPTION FOR CHANGE ]

Due to product discontinuation, the circuit protectors used in the circuit protector units (power switch) for the EM and CM have been replaced by the new ones with different appearance.

Appearance changes (both on the EM and the CM):

1. ON and OFF symbols added
2. LED color changed to white from blue
3. Switch frame changed in shape and color to black from gray

	Power OFF	Power ON
NEW		
OLD		

With this change, however, the part number and unit code of each circuit protector unit is NOT changed.

- A-EAL2-U011-00-A Circuit Protector Unit 100/115V, E:TEC6-011
- A-EAL2-U027-00-A Circuit Protector Unit 230V, E:TEC6-027
- A-EAL3-U010-00-A Circuit Protector Unit 100/115V, E:TEC6-011C
- A-EAL3-U021-00-A Circuit Protector Unit 230V, E:TEC6-021C

**[ PART CHANGE ] ( IVDP=Inventory Disposal, ICHB=Inter Changeable )**

Description	Old	QTY	IVDP	New	QTY	ICHB
N/A						

**[ REMARKS ]**

EmbeddingModule:

- J0: 5104 0312 and after
- A1: 5108 0568 and after
- E2: 5112 0429 and after
- JC2: 5116 0165 and after
- J2: 5128 0131 and after

Cryo Module:

- J0: 5105 0295 and after
- A1: 5109 0558 and after
- E2: 5113 0429 and after
- JC2: 5117 0165 and after
- J2: 5129 0131 and after

**[ DISTRIBUTION ]**

- SFA
- SFJ
- SFE
- SCN

**[ APPROVAL ]**

Name : T. Sugawara

Signature :

