



Product Support Bulletin

No. 21-01


Date: December 8, 2021

Product: Tissue-Tek® products

Subject: Addition of Alternative 20-Slide Basket

Ref. No. N/A

Please be advised that SSCL informed SFJ that they decided to use the alternative material (RENATUS) for the current 20-slide basket (4768) because of difficulty in procuring the base current material. With this change, the basket color is changed to white (below).



Left: alternative basket

Right: current basket

The quality and performance of alternative 20-slide baskets were proved equal to the current 20-slide baskets through evaluation testing by the manufacturer.

The alternative 20-slide basket is planned to be supplied from the beginning of the year 2022 as a replacement accessory item for the following instruments. In addition, the current 20-slide basket, supplied standard with these instruments, might be switched to the alternative basket sometime soon because the supply situation of the current basket material becomes worse.

Tissue-Tek Film® Automated Coverslipper

Tissue-Tek® Glas™ g2 Automated Glass Coverslipper

Tissue-Tek Prisma® Plus Automated Slide Stainer

However, please note that the change is just a temporary expedient until the situation is improved to ensure the sufficient inventory of the current basket material.

End