



# Instrument Engineering Change Notice

NO. Prisma Plus-2102

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek Prisma Plus      DATE : January 18, 2021  
 PRODUCT CODES : 6170, 6171, 6172, 6173,      ORIGINATOR : Y. Yazawa  
 IECN TYPE : DOCUMENT CHANGE      REFERENCE NUMBER : 17148  
 SERIAL NUMBERS AFFECTED : N/A      ATTACHMENTS : AK6-QE-055-05\_Prisma Plus Unpacking & Installation Manual  
 FIELD RETROFIT REQUIRED : NO

## Title: Revision of Prisma Plus Unpacking and Installation Manual (English)

### [ REASON AND DESCRIPTION FOR CHANGE ]

The Unpacking and Installation Manual was updated to include the following additions which the manufacturer had found it necessary to make in the course of updating their risk management files for QSR compliance.

- \* Set the usage management of the activated carbon filter to "In Use" and change the usage time limit as necessary.
- \* Ensure a clearance of at least 100mm between the back of the instrument and a wall.

### [ PART CHANGE ] ( IVDP=Inventory Disposal, ICHB=Inter Changeable )

Description	Old	QTY	IVDP	New	QTY	ICHB
Unpacking and Installation Manual (English)	AK6-QE-055-4	1		AK6-QE-055-05	1	

### [ REMARKS ]

### [ DISTRIBUTION ]

- SFA   
 SFJ   
 SFE   
 SCN

### [ APPROVAL ]

Name : T. Sugawara

Signature :