



Instrument Engineering Change Notice

NO. VIP6AI-2003

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek VIP6 AI DATE : July 31, 2020
PRODUCT CODES : 6040, 6041, 6042, 6043 ORIGINATOR : T. Miyasaka
IECN TYPE : PART DESIGN CHANGE REFERENCE NUMBER : 16765-1, 16800
SERIAL NUMBERS AFFECTED : See Remarks ATTACHMENTS : N/A
FIELD RETROFIT REQUIRED : NO

Title: Change of Fume Fan

[REASON AND DESCRIPTION FOR CHANGE]

The fume fan was changed to a new one due to product discontinuation.

Because the old and new fume fans are installed in a different way, the old fume fan needs to be replaced with the new one together with the fan cover as a set to maintain interchangeability. To attach the new fume fan to the new fan cover, the new and longer washer-assembled pan head screws are needed.

For convenience in replacement, the fume fan replacement unit, consisting of the new fume fan, fan cover and screws, is available as a maintenance part.

[PART CHANGE] (IVDP=Inventory Disposal, ICHB=Inter Changeable)

Description	Old	QTY	IVDP	New	QTY	ICHB
Fume Fan / E:VIP6-033	F5273200	1	NO	---	---	NO
Fume Fan / E:VIP6-053	---	---		A-EO70-U873-00-A	1	NO
Fan Cover	O70-720-00	1	NO	A-MO70-0200-00-A	1	NO
Pan Head Screw, Sems B, M4x35, SUS	B6292019	2	NO	---	---	
Pan Head Screw, Sems B, M4x40, SUS	---	---		B6292013	2	
Fume Fan Replacement Unit	---	---		A-MO70-U200-00-A	1	

[REMARKS]

Implementation:

A1: 6040 0748 and after

J0: 6041 0498 and after

E2: 6042 0757 and after

J2: 6043 0222 and after

[DISTRIBUTION]

SFA

SFJ

SFE

SCN

[APPROVAL]

Name : T. Sugawara

Signature :

