



Instrument Engineering Change Notice

NO. VIP6AI-2103

Sakura Finetek Japan Co. Ltd.

PRODUCT NAME : Tissue-Tek VIP6 AI DATE : March 5, 2021
 PRODUCT CODES : 6040, 6041, 6042, 6043 ORIGINATOR : Yazawa
 IECN TYPE : PART DESIGN CHANGE REFERENCE NUMBER : 17161-1
 SERIAL NUMBERS AFFECTED : See Remarks ATTACHMENTS : N/A
 FIELD RETROFIT REQUIRED : NO

Title: Change of Serial Number Labels (New CI)

[REASON AND DESCRIPTION FOR CHANGE]

The serial number labels were modified to uniform label layout and font type for new corporate identity. Also, the ZIP code for SFJ was changed in some labels.

For details, see the figures of the new labels on the next page.

(NOTE) This IECN only notifies of change of the serial number labels. They are not supplied as a part for service/maintenance.

[PART CHANGE] (IVDP=Inventory Disposal, ICHB=Inter Changeable)

Description	Old	QTY	IVDP	New	QTY	ICHB
Serial number label (A1)	O79-120-00	1		O79-120-01	1	YES
Serial number label (J0)	O79-121-00	1		O79-121-01	1	YES
Serial number label (E2)	O79-122-00	1		O79-122-01	1	YES
Serial number label (J2)	A-MO79-0123-01-A	1		A-MO79-0123-02-A	1	YES

[REMARKS]

Implementation:

A1: 6040 0829 and after

J0: 6041 0530 and after

E2: 6042 0859 and after

J2: 6043 0270 and after

[DISTRIBUTION]

SFA
 SFJ
 SFE
 SCN

[APPROVAL]

Name : T. Sugawara

Signature :

New Serial Number Labels

- A1 model

Tissue-Tek VIP® 6 AI
Vacuum Infiltration Processor

CODE	MODEL	VOLTS	AMPS	HERTZ
6040	VIP 6 AI-A1	115	12.0	50/60

Serial Number : 60400000-0000

Manufactured for:
Sakura Finetek USA, Inc., Torrance, CA 90501 U.S.A.
Sakura Finetek Europe B.V., 2408 AV Alphen aan den Rijn, NL.
Sakura Finetek Japan Co., Ltd., Tokyo, 103-0007, Japan

Manufactured by:
Sakura Seiki Co., Ltd., Nagano, 387-0015, Japan

O79-120-01
MADE IN JAPAN

IVD

- J0 model

ティッシュ・テック® VIP™ 6 AI 6041

型 式	VIP 6 AI-J0		
電 源 入 力	12.0A/AC100V 50/60Hz		
製造販売届出番号	20B2X00014000032	一般医療機器(クラスI)	
製 造 番 号	60410000-0000	特定保守管理医療機器	

製造販売業者 **サクラ精機株式会社** **EMC適合**
長野県千曲市大字八幡 1122-8

販売業者 **サクラファインテックジャパン株式会社** O79-121-01
東京都中央区日本橋本町 3-1-9

- E2 model

Tissue-Tek VIP® 6 AI
Vacuum Infiltration Processor

CODE	MODEL	VOLTS	AMPS	HERTZ
6042	VIP 6 AI-E2	230	6.0	50/60

Serial Number : 60420000-0000

Manufactured for:
Sakura Finetek USA, Inc., Torrance, CA 90501 U.S.A.
Sakura Finetek Europe B.V., 2408 AV Alphen aan den Rijn, NL.
Sakura Finetek Japan Co., Ltd., Tokyo, 103-0007, Japan

Manufactured by:
Sakura Seiki Co., Ltd., Nagano, 387-0015, Japan

O79-122-01
MADE IN JAPAN

IVD

CE

- J2 model

自动组织脱水机

■ 医疗器械备案凭证编号: 国械备20161905号 ■ 生产者地址:
75-5,IMOJIYA,CHIKUMA-CITY,NAGANO.PREF.,JAPAN

■ 产品名称: 自动组织脱水机

■ 产品编号: 60430000-Q000 ■ 生产场所地址:
75-5,IMOJIYA,CHIKUMA-CITY,NAGANO.PREF.,JAPAN

■ 规格型号: VIP 6 AI-J2

■ 尺寸: 606(W)×675(D)×1329(H) mm

■ 电源连接条件:
● 电流: 10A
● 电压: 230VAC ± 10%, 50/60Hz
● 湿度: 30-85% (无冷凝)
● 温度: 10-40 °C

■ 生产者名称: Sakura Seiki Co., Ltd.

■ 生产日期: 0000年00月

■ 有效期: 八年

※其他内容详见说明书

SAKURA

代理人: 樱花医疗科技(泰州)有限公司
地 址: 泰州市药城大道1号R09幢
邮编: 225300 电话: 0523-86201011 传真: 0523-86201022

M079-0123-02